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

\*SWEDAMP: Årsrapport 2022

\*SPARG: A Survey of the Lower Limb Amputee Population in Scotland, 2011

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



- **3.2 Interface components**

Interface components of a prosthesis are in direct contact with the wearer. They are the means of achieving axial stabilization, transverse stabilization and suspension (see clause 3).

Interface components include

- a) the socket,
  - b) suspensory components, which maintain the prosthesis relative to the body,
  - c) the liner.
- 
- NOTE Interface components can contain elements which contribute to the activation and/or control of functional components.

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- ISO har rekommendation om hur gränssnittet ska beskrivas (**ISO 13405-2: Prosthetics and orthotics — Classification and description of prosthetic components — Part 3: Description of lower limb prosthetic components**)
- Den fullständiga ISO-beskrivningen av Direct Socket, finns publicerad i Canadian Journal of Prosthetics and Orthotics  
<https://jps.library.utoronto.ca/index.php/cpoj/article/view/34672>

## The Transfemoral interface of the Direct Socket-TF

**Force-transmission properties:**

- The AXIAL stabilization:** is generated by a total socket stabilization conforming to the shape of the femur and soft tissues about the femur
- The TRANSVERSE stabilization:** the anteroposterior, mediolateral and rotational stabilization are modified during the direct lamination process
- SUSPENSION:** to minimize the axial movement a Seal-In liner (partial distal vacuum) is preferably used, in some cases Locking liner with Lanyard or Pin can be used
- STIFNESS:** the socket is partly flexible, i.e. the larger part of the socket is laminated with a distal 4-hole adapter and the proximal part is made of a flexible silicone. That part acts as a stabilisation for the hip joint muscles (contraction) during Initial contact, Loading responds, Mid stance and Terminal stance

# The Transfemoral interface of the Direct Socket-TF

Force-  
transmission  
properties:

The AXIAL stabilization:  
is generated by a total socket stabilization conforming to the shape of the femur and soft tissues about the femur

The TRANSVERSE stabilization:  
the anteroposterior, mediolateral and rotational stabilization are modified during the direct lamination process

SUSPENSION: to minimize the axial movement a Seal-In liner (partial distal vacuum) is preferably used, in some cases Locking liner with Lanyard or Pin can be used

STIFNESS: the socket is partly flexible, i.e. the larger part of the socket is laminated with a distal 4-hole adapter and the proximal part is made of a flexible silicone. That part acts as a stabilisation for the hip joint muscles (contraction) during Initial contact, Loading responds, Mid stance and Terminal stance

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

# Forskningsartiklar

**CPDJ** OPEN ACCESS  
Journal Homepage: <https://doi.org/10.33137/cpj.0101>  
Volume 3, Issue 2, Article No. 6, 2020

**RESEARCH ARTICLE**

**TRANSFEMORAL SOCKET FABRICATION METHOD USING DIRECT CASTING: OUTCOMES REGARDING PATIENT SATISFACTION WITH DEVICE AND SERVICES**

Marable W.R.<sup>1</sup>, Smith C.<sup>1</sup>, Sigurjónsson B.P.<sup>2</sup>, Atlason I.F.<sup>3</sup>, Johannesson G.A.<sup>2,4\*</sup>

<sup>1</sup> Osseur Americas, Foothill Ranch, California, USA.  
<sup>2</sup> Osseur HF, Reykjavik, Iceland.  
<sup>3</sup> Quik Lookse, Reykjavik, Iceland.  
<sup>4</sup> TeamÖlmed, Stockholm, Sweden.

**ABSTRACT**

**BACKGROUND:** Direct Socket for transfemoral (DS-TF) prosthetic user is a novel method of fabricating a laminated interface on to the residual limb but requires different training, production method and service model than what most prosthetists are familiar with. This method and model may improve patient satisfaction by enabling interface fabrication and delivery in one visit.

**OBJECTIVES:** Document patient satisfaction regarding DS-TF interface versus the prosthetic users' previous socket in terms of interface function and the clinic service model.

**METHODOLOGY:** In this longitudinal study (from July 2018 to April 2020), the DS-TF was implemented in six prosthetic clinics across the United States. Certified prosthetists (CP) and assistants were trained using a standard protocol. 47 prosthetic users participated, both those in need of a new socket and those without need. Two modules from the Orthotics and Prosthetics Users' Survey (OPUS), involving questions related to satisfaction with the Device and Services, was used to evaluate each DS-TF user outcome vs. baseline. The only part of the prosthesis that was replaced was the interface, except in 2 cases.

**FINDINGS:** Each DS-TF interface was fabricated, fit and delivered in a single clinic visit. At 6-month follow-up, 38 users reported an average of 20.8% increase in satisfaction with their new interface compared with original, and a 14.8% increase in satisfaction with the services they received from the clinic in providing of the new prosthesis vs. their original prosthesis. The main outcome increases were between baseline (initial fitting) and 6-week follow-up and remained consistent after 6 months. This improvement was consistent irrespective if the user needed a new socket for clinical reasons or not.

**CONCLUSION:** This study shows that after a standardized training and implementation, the DS-TF fabrication process including a new interface improves the user's satisfaction with their prosthetic device and services.

**INTRODUCTION**

Transfemoral (TF) amputation(s) is a devastating procedure for any person and expensive in terms of acute healthcare cost<sup>1</sup> and rehabilitation cost.<sup>1,2</sup> TF amputation highly restricts the amputee's overall mobility<sup>3</sup> and requires a good functional prosthesis, especially the interface (i.e. the socket, the liner and its suspension) to enable the amputee to regain as much of their previous mobility as possible.<sup>4,5</sup> TF amputation is also associated with longer amputee rehabilitation time compared with transtibial (TT) amputation and more need of assistance, especially in elderly patients.<sup>1</sup> TF amputees report their general health related quality of life to be lower than that of non-amputees or TT amputees. Furthermore, specific problems have been related to the use of TF prostheses, including limited hip joint range of motion that restricts users comfort or ability to perform daily activities such as sitting, walking, picking

**ARTICLE INFO**

Received: August 3, 2020  
Accepted: November 16, 2020  
Published: November 23, 2020

**CITATION**

Marable W.R., Smith C., Sigurjónsson B.P., Atlason I.F., Johannesson G.A. Transfemoral socket fabrication method using direct casting: outcomes regarding patient satisfaction with device and services. Canadian Prosthetics & Orthotics Journal. 2020;Volume 3, Issue 2, No.6. <https://doi.org/10.33137/cpj.v3i2.34872>

**KEYWORDS**

Transfemoral amputation, Amputee, Prostheses, Socket, Interface, Outcome measure, Satisfaction, Direct casting

\* CORRESPONDING AUTHOR:  
G. Anton Johannesson, PhD  
TeamÖlmed, Kättinggatan 12, 164 40 Kista, Stockholm, Sweden.  
E-mail: [gjohannesson@teamolmed.se](mailto:gjohannesson@teamolmed.se)  
ORCID: <https://orcid.org/0000-0001-8726-450X>

Marable W.R., Smith C., Sigurjónsson B.P., Atlason I.F., Johannesson G.A. Transfemoral socket fabrication method using direct casting: outcomes regarding patient satisfaction with device and services. Canadian Prosthetics & Orthotics Journal. 2020;Volume 3, Issue 2, No.6. <https://doi.org/10.33137/cpj.v3i2.34872>

- Users' satisfaction OPUS
  - Device
  - Services

- Quality of life EQ-5D-5L
- Mobility Plus-M
- Stability/Suspension CLASS
- Balance confidence ABC

Marable W.R, Smith C, Sigurjónsson B.P, Atlason I.F, Johannesson G.A. Transfemoral socket fabrication method using direct casting: outcomes regarding patient satisfaction with device and services. Canadian Prosthetics & Orthotics Journal. 2020;Volume 3, Issue 2, No.6.

**CPDJ** OPEN ACCESS  
Journal Homepage: <https://doi.org/10.33137/cpj.0101>  
Volume 4, Issue 1, Article No. 6, 2021

**RESEARCH ARTICLE**

**CLINICAL OUTCOME OF TRANSFEMORAL DIRECT SOCKET INTERFACE (PART 2)**

Walker J.<sup>1</sup>, Marable W.R.<sup>2,3</sup>, Smith C.<sup>1</sup>, Sigurjónsson B.P.<sup>2</sup>, Atlason I.F.<sup>1</sup>, Johannesson G.A.<sup>2,4\*</sup>

<sup>1</sup> Virginia Prosthetic & Orthotic, Roanoke, Virginia, USA.  
<sup>2</sup> Osseur HF, Foothill Ranch, California, USA.  
<sup>3</sup> Osseur HF, Reykjavik, Iceland.  
<sup>4</sup> Quik Lookse, Reykjavik, Iceland.  
<sup>5</sup> TeamÖlmed, Stockholm, Sweden.

**ABSTRACT**

**BACKGROUND:** Amputation at the transfemoral (TF) level reduces the rate of successful prosthetic fitting, functional outcome, and quality of life (QoL) compared with transtibial amputation. The TF socket interface is considered the most critical part of the prosthesis, but socket discomfort is still the most common user complaint. Direct Socket for transfemoral prosthesis users is a novel interface fabrication process where the socket is shaped and laminated directly on the residual limb and delivered in a single visit.

**OBJECTIVES:** The aim of this study was to investigate if prosthetic users' quality of life (QoL), comfort, and mobility with a Direct Socket TF interface were comparable to their experience with their previous prostheses.

**METHODOLOGY:** The pre/post design prospective cohort study included 47 subjects. From this cohort, 36 subjects completed the 6-month follow-up (mean age 58 years, 27 males). Outcomes at baseline included EQ-5D-5L\*, PLUS-M\*\*, CLASS, ABC, AMPPRO, and TUQ. At 6-weeks and 6-months, subjects repeated all measures. Seven Certified Prosthetist (CP) investigators performed observations and data collection at six different sites (from July 2018 to April 2020).

**FINDINGS:** Results showed significant improvement in all outcome measures for the 36 subjects that completed both 6-weeks and 6-month follow-ups. CLASS sub-scales showed significantly improved stability, suspension, comfort, and socket appearance. Improvement in K-Level and less use of assistive devices were observed with the AMPPRO instrument, indicating improved user mobility and performance. QoL was also increased, as measured in Quality-Adjusted-Life-Years (QALY) from the EQ-5D-5L.

**CONCLUSION:** Evidence from the findings demonstrate that the Direct Socket TF system and procedure can be a good alternative to the traditional method of prosthetic interface delivery.

**ARTICLE INFO**

Received: February 19, 2021  
Accepted: May 26, 2021  
Published: June 8, 2021

**CITATION**

Walker J, Marable W.R, Smith C, Sigurjónsson B.P, Atlason I.F, Johannesson G.A. Clinical outcome of transfemoral direct socket interface (part 2). Canadian Prosthetics & Orthotics Journal. 2021;Volume 4, Issue 1, No.6. <https://doi.org/10.33137/cpj.v4i1.30005>

**KEYWORDS**

Transfemoral Amputation, Health, Amputee, Prostheses, Socket, Interface, Comfort, Outcome Measure, Satisfaction, Direct Casting

**INTRODUCTION**

The primary goal for people that have undergone lower limb amputation is to return to the main activities of daily living, including recreational and professional activities.<sup>1</sup> Transfemoral (TF) amputation and the subsequent loss of knee function is known to negatively impact prosthetic fitting, functional outcome, and quality of life (QoL), as compared to transtibial amputation.<sup>2,3</sup> Persons with TF amputation face significantly more challenges when receiving a prosthesis for the first time. With more days spent in rehabilitation, the functional outcome and QoL are still lower at the time of discharge compared to persons with transtibial (TT) amputation.<sup>3,4</sup> Additionally, the Certified Prosthetist (CP) faces a more significant challenge in fitting TF patients versus TT patients.<sup>3,13</sup> TF amputations may account for approximately 40% of all lower-limb amputations in the US alone and result in a twofold higher mortality rate than after TT amputations<sup>15</sup> and a lower rate

\* CORRESPONDING AUTHOR:  
G. Anton Johannesson, PhD  
TeamÖlmed, Kättinggatan 12, 164 40 Kista, Stockholm, Sweden.  
E-mail: [gjohannesson@teamolmed.se](mailto:gjohannesson@teamolmed.se)  
ORCID: <https://orcid.org/0000-0001-8726-450X>

Walker J, Marable W.R, Smith C, Sigurjónsson B.P, Atlason I.F, Johannesson G.A. Clinical outcome of transfemoral direct socket interface (part 2). Canadian Prosthetics & Orthotics Journal. 2021;Volume 4, Issue 1, No.6. <https://doi.org/10.33137/cpj.v4i1.30005>

Walker J, Marable W.R, Smith C, Sigurjónsson B.P, Atlason I.F, Johannesson G.A. clinical outcome of transfemoral direct socket interface (part 2). Canadian Prosthetics & Orthotics Journal. 2021;Volume 4, Issue 1, No.6.

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29.8% INCREASE IN  
SATISFACTION WITH  
INTERFACE AFTER  
6 MONTHS

14.8% INCREASE  
IN SATISFACTION  
WITH SERVICE

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